



TECHNICAL DATA


<i>Voltage</i>	<i>nominal</i>	125	<i>kV</i>
	<i>maximum for test</i>	140	
<i>Inverse voltage</i>	<i>nominal</i>	125	<i>kV</i>
	<i>maximum for test</i>	140	
<i>Focal spot (IEC 60336:2005)</i>		0.6	<i>mm</i>
<i>Filament characteristics</i>		3.0 ÷ 6.0	<i>V</i>
		3.0 ÷ 4.0	<i>A</i>
<i>Filament frequency limit</i>		20	<i>kHz</i>
<i>Anode material</i>		Tungsten	
<i>Target angle</i>		9	<i>°</i>
<i>Anode heat storage capacity</i>		60000	<i>J</i>
<i>Maximum anode cooling rate</i>		600	<i>W</i>
<i>Nominal anode input power at 1.0 s (DC)</i>		2020	<i>W</i>
<i>Inherent filtration</i>		0.5 mm Al (±0.05 at 75 kV)	
<i>Maximum diameter</i>		49	<i>mm</i>
<i>Overall length</i>		173	<i>mm</i>

WORKING CONDITIONS


<i>Approx Weight</i>	660 g (1230 g with lead shield)
<i>X-ray coverage</i>	Ø285 mm at 1000 mm SID
<i>Cooling method</i>	Oil immersed (60°C max) convection cooling
<i>Tube holding</i>	Holding the screw of the anode shank, or holding the glass at anode and cathode side with insulating material
<i>Operating limits (dielectric oil)</i>	
<i>Oil temperature</i>	15 ÷ 60 °C
<i>Oil pressure</i>	70 ÷ 140 kPa
<i>Shipping and storage limits</i>	
<i>Temperature</i>	-40 ÷ 70 °C
<i>Humidity (no condensation)</i>	10 ÷ 90 %
<i>Atmospheric pressure</i>	50 ÷ 110 kPa
<i>Seasoning / Schedule</i>	
The tube needs seasoning if not used for 3 months or at the first usage. For conditioning procedure start at the lowest kV, at max 3 mA, exposures of 30 sec. Rise kV till maximum at step of maximum 10 kV. If arch happens restart from previous kV step.	

DIMENSIONS

STANDARD



LEAD SHIELD



GENERAL INFORMATIONS

When mounting tube inserts adopt proper caution, in order to avoid glass bulb breaking and fragments projection. Please use protective gloves and glasses.
 Tube insert connected to H.V. supply is a radiation source: be sure to take all necessary safety cautions


- Wash thoroughly with alcohol the external surface of tube insert (care of fire risk). Avoid contact of dirty surfaces with cleaned tube insert.
- Clamp system inside housing or self-contained units must not mechanically stress the tube.
- After installation, check the right working of the tube (no fluctuation of tube current nor crackling)
- Comply with insert thermal parameters, planning and programming the exposure parameters and cooling pauses. Housing or self-contained units must be provided with an adequate thermic protection.
- Voltages indicated in charts are valid for transformer supplied with ground center.
- Tube inserts contain environment polluting materials, particularly lead liner tubes. Please apply to qualified operator for waste disposal, according to local regulation requirements.

INCIDENT REPORT ACCORDING TO 93/42/EEC MEDICAL DEVICES DIRECTIVE


In order to comply with CE marking requirements, end users have to report to local Competent Authority all the informations about possible incidents involving the device, regarding any deterioration in its characteristics and performances, as well as any inaccuracies in this documentation, which might lead to or might have led to the death of patient / user or a deterioration in his state of health. This information must be promptly reported also to C.E.I. in order to start manufacturer reporting, as per above mentioned directive.




EMISSION CURVES DC




FILAMENT CHARACTERISTICS



THERMAL CURVES




RATING CHARTS DC





RATING CHARTS DC 25 fps - 50% DUTY CYCLE



RATING CHARTS DC 12.5 fps - 50% DUTY CYCLE

